

Anatomic Healing Abutment™

by MEGA'GEN

Optimal Gingiva Healing

Aesthetic final prosthesis,
like natural teeth



ANYRIDGE®
by MEGA'GEN

Anatomic Healing Abutments



Healing Abutment

(Anatomic type)

- Use with a Hand Driver(1.2 Hex).
- Abutment Screw included.H=4 ARHAS1804/
H=5 ARHAS1805/ H=7 ARHAS1807
- Used for non-submerged type surgery or for two stage surgery.
- Choose appropriate diameter and height of Healing Abutment according to situation.
- Helps to form suitable emergence profile during period of gingival healing.
- Recommend torque : by hand (5 - 8Ncm)

Type	MD (mm)	LL (mm)	Height (mm)	Connection	Ref.C
Incisor	4.0	5.0	4	Hex	ARHI40504T
			5		ARHI40505T
			7		ARHI40507T
	4.5	4.5	4		ARHI45454T
			5		ARHI45455T
			7		ARHI45457T
	6.0	5.0	4		ARHI60504T
			5		ARHI60505T
			7		ARHI60507T
	7.0	6.0	4	Non-Hex	ARHI70604T
			5		ARHI70605T
			7		ARHI70607T
	4.0	5.0	4		ARHI40504NT
			5		ARHI40505NT
			7		ARHI40507NT
	4.5	4.5	4		ARHI45454NT
			5		ARHI45455NT
			7		ARHI45457NT
	6.0	5.0	4		ARHI60504NT
			5		ARHI60505NT
			7		ARHI60507NT
	7.0	6.0	4		ARHI70604NT
			5		ARHI70605NT
			7		ARHI70607NT



Type	MD (mm)	LB (mm)	Height (mm)	Connection	Ref.C
Canine	5.0	6.5	4	Hex	ARHC50654T
			5		ARHC50655T
			7		ARHC50657T
	5.0	6.5	4	Non-Hex	ARHC50654NT
			5		ARHC50655NT
			7		ARHC50657NT



Type	MD (mm)	LB (mm)	Height (mm)	Connection	Ref.C
Pre-Molar	4.5	6.0	4	Hex	ARHM45604T
			5		ARHM45605T
			7		ARHM45607T
	5.0	7.0	4		ARHM50704T
			5		ARHM50705T
			7		ARHM50707T
	4.5	6.0	4	Non-Hex	ARHM45604NT
			5		ARHM45605NT
			7		ARHM45607NT
	5.0	7.0	4		ARHM50704NT
			5		ARHM50705NT
			7		ARHM50707NT



Type	MD (mm)	LB (mm)	Height (mm)	Connection	Ref.C
Molar	6.0	7.0	4	Hex	ARHM60704T
			5		ARHM60705T
			7		ARHM60707T
	6.0	8.0	4		ARHM60804T
			5		ARHM60805T
			7		ARHM60807T
	6.0	9.0	4		ARHM60904T
			5		ARHM60905T
			7		ARHM60907T
	7.0	8.0	4		ARHM70804T
			5		ARHM70805T
			7		ARHM70807T
	7.0	9.0	4		ARHM70904T
			5		ARHM70905T
			7		ARHM70907T
	7.0	10.0	4		ARHM70104T
			5		ARHM70105T
			7		ARHM70107T
	8.0	9.0	4		ARHM80904T
			5		ARHM80905T
			7		ARHM80907T
	8.0	10.0	4		ARHM80104T
			5		ARHM80105T
			7		ARHM80107T
	6.0	7.0	4	Non-Hex	ARHM60704NT
			5		ARHM60705NT
			7		ARHM60707NT
	6.0	8.0	4		ARHM60804NT
			5		ARHM60805NT
			7		ARHM60807NT
	6.0	9.0	4		ARHM60904NT
			5		ARHM60905NT
			7		ARHM60907NT
	7.0	8.0	4		ARHM70804NT
			5		ARHM70805NT
			7		ARHM70807NT
	7.0	9.0	4		ARHM70904NT
			5		ARHM70905NT
			7		ARHM70907NT
	7.0	10.0	4		ARHM70104NT
			5		ARHM70105NT
			7		ARHM70107NT
	8.0	9.0	4		ARHM80904NT
			5		ARHM80905NT
			7		ARHM80907NT
	8.0	10.0	4		ARHM80104NT
			5		ARHM80105NT
			7		ARHM80107NT



Type	MD (mm)	LB (mm)	Height (mm)	Connection	Ref.C
Special	4.5	6.0	4	Hex	ARHS45604T
			5		ARHS45605T
			7		ARHS45607T
	5.0	6.5	4		ARHS50654T
			5		ARHS50655T
			7		ARHS50657T
	5.0	7.0	4		ARHS50704T
			5		ARHS50705T
			7		ARHS50707T
	6.0	7.0	4		ARHS60704T
			5		ARHS60705T
			7		ARHS60707T
	6.0	8.0	4		ARHS60804T
			5		ARHS60805T
			7		ARHS60807T
	6.0	9.0	4		ARHS60904T
			5		ARHS60905T
			7		ARHS60907T
	7.0	8.0	4		ARHS70804T
			5		ARHS70805T
			7		ARHS70807T
	7.0	9.0	4		ARHS70904T
			5		ARHS70905T
			7		ARHS70907T
	7.0	10.0	4		ARHS70104T
			5		ARHS70105T
			7		ARHS70107T
	8.0	9.0	4		ARHS80904T
			5		ARHS80905T
			7		ARHS80907T
	8.0	10.0	4		ARHS80104T
			5		ARHS80105T
			7		ARHS80107T
	4.5	6.0	4	Non-Hex	ARHS45604NT
			5		ARHS45605NT
			7		ARHS45607NT
	5.0	6.5	4		ARHS50654NT
			5		ARHS50655NT
			7		ARHS50657NT
	5.0	7.0	4		ARHS50704NT
			5		ARHS50705NT
			7		ARHS50707NT
	6.0	7.0	4		ARHS60704NT
			5		ARHS60705NT
			7		ARHS60707NT
	6.0	8.0	4		ARHS60804NT
			5		ARHS60805NT
			7		ARHS60807NT
	6.0	9.0	4		ARHS60904NT
			5		ARHS60905NT
			7		ARHS60907NT
	7.0	8.0	4		ARHS70804NT
			5		ARHS70805NT
			7		ARHS70807NT
	7.0	9.0	4		ARHS70904NT
			5		ARHS70905NT
			7		ARHS70907NT
	7.0	10.0	4		ARHS70104NT
			5		ARHS70105NT
			7		ARHS70107NT
	8.0	9.0	4		ARHS80904NT
			5		ARHS80905NT
			7		ARHS80907NT
	8.0	10.0	4		ARHS80104NT
			5		ARHS80105NT
			7		ARHS80107NT

BLUEDIAMOND IMPLANT

Anatomic Healing Abutments

Healing Abutment

(Anatomic type)

- Abutment screw (H=4 AROHAS1604/ H=5 AROHAS1605/ H=7 AROHAS1607) included

- For use with non-submerged type surgery & 2-stage surgery
- Different diameter & height options
- Helps to form suitable emergence profile during gingival healing
- Recommend torque: by hand (5 - 8Nm)
- Use with Hand Driver(1.2 Hex)




Type	MD (mm)	LL (mm)	Height (mm)	Connection	Ref.C
Incisor	4.0	5.0	4	Octa	AROHIN40504T
			5		AROHIN40505T
			7		AROHIN40507T
	4.5	4.5	4		AROHIN45454T
			5		AROHIN45455T
			7		AROHIN45457T
	6.0	5.0	4		AROHIN60504T
			5		AROHIN60505T
			7		AROHIN60507T
	7.0	6.0	4		AROHIN70604T
			5		AROHIN70605T
			7		AROHIN70607T
	4.0	5.0	4	Non-Octa	AROHIN40504NT
			5		AROHIN40505NT
			7		AROHIN40507NT
	4.5	4.5	4		AROHIN45454NT
			5		AROHIN45455NT
			7		AROHIN45457NT
	6.0	5.0	4		AROHIN60504NT
			5		AROHIN60505NT
			7		AROHIN60507NT
	7.0	6.0	4		AROHIN70604NT
			5		AROHIN70605NT
			7		AROHIN70607NT

Type	MD (mm)	LL (mm)	Height (mm)	Connection	Ref.C
Incisor	4.0	5.0	4	Octa	AROHIR40504
			5		AROHIR40505
			7		AROHIR40507
	4.5	4.5	4		AROHIR45454
			5		AROHIR45455
			7		AROHIR45457
	6.0	5.0	4		AROHIR60504
			5		AROHIR60505
			7		AROHIR60507
	7.0	6.0	4		AROHIR70604
			5		AROHIR70605
			7		AROHIR70607
	4.0	5.0	4	Non-Octa	AROHIR40504N
			5		AROHIR40505N
			7		AROHIR40507N
	4.5	4.5	4		AROHIR45454N
			5		AROHIR45455N
			7		AROHIR45457N
	6.0	5.0	4		AROHIR60504N
			5		AROHIR60505N
			7		AROHIR60507N
	7.0	6.0	4		AROHIR70604N
			5		AROHIR70605N
			7		AROHIR70607N



Type	MD (mm)	LB (mm)	Height (mm)	Connection	Ref.C
Canine	5.0	5.5	4	Octa	AROHCN50654T
			5		AROHCN50655T
			7		AROHCN50657T
	5.0	5.5	4	Non-Octa	AROHCN50654NT
			5		AROHCN50655NT
			7		AROHCN50657NT

Type	MD (mm)	LB (mm)	Height (mm)	Connection	Ref.C
Canine	5.0	5.5	4	Octa	AROHC50654T
			5		AROHC50655T
			7		AROHC50657T
	5.0	5.5	4	Non-Octa	AROHC50654NT
			5		AROHC50655NT
			7		AROHC50657NT



Type	MD (mm)	LB (mm)	Height (mm)	Connection	Ref.C
Pre-Molar	4.5	6.0	4	Octa	AROHMR45604T
			5		AROHMR45605T
			7		AROHMR45607T
	5.0	7.0	4		AROHMR50704T
			5		AROHMR50705T
			7		AROHMR50707T
	4.5	6.0	4	Non-Octa	AROHMR45604NT
			5		AROHMR45605NT
			7		AROHMR45607NT
	5.0	7.0	4		AROHMR50704NT
			5		AROHMR50705NT
			7		AROHMR50707NT



Type	MD (mm)	LB (mm)	Height (mm)	Connection	Ref.C
Special	4.5	6.0	4	Octa	AROHSR45604T
			5		AROHSR45605T
			7		AROHSR45607T
	5.0	6.5	4		AROHSR50654T
			5		AROHSR50655T
			7		AROHSR50657T
	5.0	7.0	4		AROHSR50704T
			5		AROHSR50705T
			7		AROHSR50707T
	6.0	7.0	4		AROHSR60704T
			5		AROHSR60705T
			7		AROHSR60707T
	6.0	8.0	4		AROHSR60804T
			5		AROHSR60805T
			7		AROHSR60807T
	6.0	9.0	4		AROHSR60904T
			5		AROHSR60905T
			7		AROHSR60907T
	7.0	8.0	4		AROHSR70804T
			5		AROHSR70805T
			7		AROHSR70807T
	7.0	9.0	4		AROHSR70904T
			5		AROHSR70905T
			7		AROHSR70907T
	7.0	10.0	4		AROHSR70104T
			5		AROHSR70105T
			7		AROHSR70107T
Special	8.0	9.0	4	Non-Octa	AROHSR80904T
			5		AROHSR80905T
			7		AROHSR80907T
	8.0	10.0	4		AROHSR80104T
			5		AROHSR80105T
			7		AROHSR80107T
	4.5	6.0	4		AROHMR45604NT
			5		AROHMR45605NT
			7		AROHMR45607NT
	5.0	6.5	4		AROHMR50654NT
			5		AROHMR50655NT
			7		AROHMR50657NT
	5.0	7.0	4		AROHMR50704NT
			5		AROHMR50705NT
			7		AROHMR50707NT
	6.0	7.0	4		AROHMR60704NT
			5		AROHMR60705NT
			7		AROHMR60707NT
	6.0	8.0	4		AROHMR60804NT
			5		AROHMR60805NT
			7		AROHMR60807NT
	6.0	9.0	4		AROHMR60904NT
			5		AROHMR60905NT
			7		AROHMR60907NT
	7.0	8.0	4		AROHMR70804NT
			5		AROHMR70805NT
			7		AROHMR70807NT
	7.0	9.0	4		AROHMR70904NT
			5		AROHMR70905NT
			7		AROHMR70907NT
	7.0	10.0	4		AROHMR70104NT
			5		AROHMR70105NT
			7		AROHMR70107NT
	8.0	9.0	4		AROHMR80904NT
			5		AROHMR80905NT
			7		AROHMR80907NT
	8.0	10.0	4		AROHMR80104NT
			5		AROHMR80105NT
			7		AROHMR80107NT

Type	MD (mm)	LB (mm)	Height (mm)	Connection	Ref.C
Molar	6.0	7.0	4	Octa	AROHMR60704T
			5		AROHMR60705T
			7		AROHMR60707T
	6.0	8.0	4		AROHMR60804T
			5		AROHMR60805T
			7		AROHMR60807T
	6.0	9.0	4		AROHMR60904T
			5		AROHMR60905T
			7		AROHMR60907T
	7.0	8.0	4		AROHMR70804T
			5		AROHMR70805T
			7		AROHMR70807T
	7.0	9.0	4		AROHMR70904T
			5		AROHMR70905T
			7		AROHMR70907T
	7.0	10.0	4		AROHMR70104T
			5		AROHMR70105T
			7		AROHMR70107T
	8.0	9.0	4		AROHMR80904T
			5		AROHMR80905T
			7		AROHMR80907T
	8.0	10.0	4		AROHMR80104T
			5		AROHMR80105T
			7		AROHMR80107T
	6.0	7.0	4	Non-Octa	AROHMR60704NT
			5		AROHMR60705NT
			7		AROHMR60707NT
	6.0	8.0	4		AROHMR60804NT
			5		AROHMR60805NT
			7		AROHMR60807NT
	6.0	9.0	4		AROHMR60904NT
			5		AROHMR60905NT
			7		AROHMR60907NT
	7.0	8.0	4		AROHMR70804NT
			5		AROHMR70805NT
			7		AROHMR70807NT
	7.0	9.0	4		AROHMR70904NT
			5		AROHMR70905NT
			7		AROHMR70907NT
	7.0	10.0	4		AROHMR70104NT
			5		AROHMR70105NT
			7		AROHMR70107NT
	8.0	9.0	4		AROHMR80904NT
			5		AROHMR80905NT
			7		AROHMR80907NT
	8.0	10.0	4		AROHMR80104NT
			5		AROHMR80105NT
			7		AROHMR80107NT

Anatomic Healing Abutments

Location Guide

